

Product Family

Honeywell Titmus Software

Honeywell Titmus Software is a specialised Vision Screener management software developed for the Honeywell Titmus V4 Vision Screener. A demo version of the software is available for 10 days **Download Software Version 2015.01** and **Quick Start Guide**

Honeywell

Product Numbers & Ordering Information

Product Numbers	Details
19386	TITMUS SOFTWARE & CABLE,V Software for Titmus VSeries Vision Screeners, model Titmus V4 Includes cables to connect to vision screener

Additional Information

Software Order Process:

- 1. Customer places order with our distributor. Distributor collects information from customer.
- 2. Distributor completes the **Software Order Form** and places order with us along with the Software Order Form.
- 3. We drop-ship the following to the customer directly
- Cables that connect the PC (software) to the vision screening instrument.
- Page that shows the webpage to download the software AND activation code.
- 20 days to use the activation code

Overview

Key Features

Subjects

• Manage your facilities subjects (employees) and keep track of their key Occupation Health details.

Vision

Control your vision screener directly from the PC and automatically record your subject's vision test results. Supports
Honeywell TITMUSi400 and the new Titmus V4.

Integration

• Open interfaces provide accessible integration points for subject and test information to & from other systems; e.g. SAP.

Business Intelligence

 Advanced statistics and reporting infrastructure allows you to extract all the information needed for printing, creating .pdf for e-mail, or creating .csv for manipulating the data in Excel.

Standards

• Can customize per your requirements or use the generic standards

Pass/Fail Viewing

• Can see the pass/fail results for each test on the screen, as you perform the test

Corrected / uncorrected vision testing

• Can test a person with and without glasses and compare the results

Honeywell Titmus Software

Reverse

• Can go back to an individual test within a complete test, if you want to redo the test. Example: You have completed the acuity tests and are now doing the color vision test. You can easily go back to do the acuity test again.

DEMO VERSION OF SOFTWARE

The **Demo Version**of the software is almost a fully functional version of our software. This is to help our customers review the software before purchase. The only difference between the demo version and the fully functional version is as follows –

- No cable provided with the demo version. Hence you cannot connect a vision screening instrument to the demo version.
- Section of the software that shows the instrument status is unavailable in the demo version, since it is not connected to the vision screening instrument.
- Can set up patients, do a full test, create a test report. But you cannot print a report from the demo version of the software.
- The demo version does NOT need an activation key to try out the software. You can ignore this screen and test the software.

Specifications

Memory Requirements

• Minimum 512 MB, 2GB required for large sites (more than 100,000 patients).

Operating System

Windows XP or later

Java

Oracle Sun Java Runtime 6

Can be used with Occupational Health Slide Sets only

• (Landolt, Snellen, TNO)

Software Version

• 2015.01

Literature & Documents

VS Software Order Form_rev M_2015-03

http://www.honeywellsafety.com/Supplementary/Documents_and_Downloads/Eye_and_Face_Protection /4294989809/1033.aspx

SoftwareQuickStartGuide

http://www.honeywellsafety.com/Supplementary/Documents_and_Downloads/Eye_and_Face_Protection /Vision_Screening/4294971115/1033.aspx

Additional Information

Download Honeywell Titmus Software Version 2015.01 Download Quick Start Guide

Software Version History

Technical Support

Sperian Protection Optical Inc, a Honeywell company

690 HP Way

Chester, VA 23836

Tel: 1-800-446-1802 (US) or+1-804-452-5200 (Intl)

Email: info@titmus.com

Australia/New Zealand

Email: hlshspamarketing@Honeywell.com

Honeywell Titmus Software

Ph Aus: 1300 139 166 Ph NZ: 0800 322 200

RELATED PRODUCTS

© Honeywell International Inc